



Sarah Tanja <stanja@uw.edu>

Method or protocol written out for how to measure chlorophyll with the NanoDrop 2000c

2 messages

CAD - Nanodrop Technical Support

<nanodrop@thermofisher.com>

To: "stanja@uw.edu" <stanja@uw.edu>

Tue, May 28, 2024 at
9:45 AM

Hello Sarah,

Thank you for your email. Please see link to below to a protocol that we have published for measuring chlorophyll on the NanoDrop One. It is important to note that the pre-defined custom method mentioned in the document is only compatible with the NanoDrop One, and not with the NanoDrop 2000.

[quantify-chlorophyll-a-and-chlorophyll-b-with-custom-method-T141.pdf](#)
(thermofisher.com)

You can replicate the Custom Method on your NanoDrop 2000 by following the attached protocol and using the below guidelines in the "Method Editor" application:

- Wavelength range Visible (350–850 nm)
- Extinction coefficient 74.8 g/L
- Analysis wavelength 666 nm
- Correction for analysis wavelength 750 nm
- Baseline correction 750 nm
- Automated pathlength On

If you have any questions, please do not hesitate to contact me; it is my pleasure to assist you. I can best be reached by email at nanodrop@thermofisher.com, and typing Attention: Moira in the subject bar. If I am unavailable, any member of our Technical Support Team would be happy to assist you at 877-724-7690.

Best,
Moira Abbott
Technical Application Specialist II
Thermo Scientific NanoDrop Products

Thermo Fisher Scientific
[3411 Silverside Road](#) | Tatnall Building, Suite 100 | Wilmington, DE 19810
Tel: +1 (302) 479-7707 | Toll Free: +1 (877) 724-7690 ext. 2400109 | Fax: +1 (302) 792-7155

Description

*****Arrived in ELOQUA on 5/23/2024*****

Hi!

Do you have a method or protocol written out for how to measure chlorophyll with the NanoDrop 2000c?

ref:!00D6007AQI.!5004z01woc5n:ref

Sarah Tanja <stanja@uw.edu>

Tue, May 28, 2024 at 10:21 AM

To: CAD - Nanodrop Technical Support <nanodrop@thermofisher.com>

Thank you so much for the info!

[Quoted text hidden]